

X-gal Staining of Cells

1. Fix cells in ice cold X-gal fixative solution for 4 minutes in a chemical fume hood.
The fixing itself can be at room temperature.
2. Rinse with PBS, let sit for 1 minute.
3. Rinse with more PBS, let sit for 10 minutes.
4. Dilute 40 mg/ml X-gal 1:40 into staining solution immediately before use. Replace PBS with diluted X-gal solution.
5. Incubate in a humidified chamber at 37°C (cell culture incubator is fine) for one-several hours until color has developed to a maximum.
6. Remove X-gal solution and rinse in PBS for several minutes.

Phase optics are not good for viewing X-gal stained cells because it masks the difference between lightly stained and unstained cells. A convenient way to see both stained and unstained cells is to push the phase slider part of the way in, which gives a fake Nomarski-style image. True Nomarski (D.I.C.) is also a good option if available, but if cells were grown in tissue culture plastic dishes, the plastic can interfere with the quality of the Nomarski optics.

X-gal fixative solution

4% formaldehyde

0.5% glutaraldehyde

0.1 M Na phosphate buffer, pH 7.2

Formaldehyde can be diluted from concentrated ampules (cleanest and safest source).

Alternatively, heat 200 ml of Na phosphate buffer with 16 g paraformaldehyde powder in a chemical fume hood with stirring until dissolved. Add 8 ml of 25% glutaraldehyde and Na phosphate buffer up to 400 ml. Store at 4°C for at least several months. *Note: Using formalin (37% liquid formaldehyde) instead of ampules or paraformaldehyde powder as a source of formaldehyde appears to give identical results.*

X-gal staining solution

0.82 g $\text{K}_3\text{Fe}(\text{CN})_6$ (5 mM final)

1.05 g $\text{K}_4\text{Fe}(\text{CN})_6 \cdot 3\text{H}_2\text{O}$ (5 mM final)

1 ml of 1 M MgCl_2 (2 mM final)

PBS up to 500 ml

Store in the dark at room temperature for several months or until considerable precipitate forms.

Some precipitate at the bottom is normal; do not shake the solution.